





REBOOTS PROGRAM SEQUENCES



REBOOTS ONE PRO

PROGRAM OVERVIEW

A Flow**B** Intense**C** Double**D** Impulse

"Flow" exerts a punctual pressure and is therefore good for loosening the muscles. The muscles are gradually compressed from the feet towards the body. We recommend this program especially for warm-up.

"Intense" most intensively promotes lactate removal. The sustained pressure in the lower chambers eliminates reflux. As a beginner, you should approach this program slowly. At high pressure, program B is very intense and effective!

Similar to program A, but two chambers are filled directly per step. Thus, "Double" achieves twice the number of passes in the same time and is suitable for shorter warm-up massages.

"Impulse" is for short and intense sessions. Since all chambers are pulsed at the same time, "Impulse" maximizes the number of possible compression passes in a given time.

PROGRAM OVERVIEW


E Intense Double



F Intense Flow



G Intense x Double X Impulse




"Intense Double" offers a similar sequence as program E, but two compression waves run in parallel. This enables us to provide a double-intensity application specifically for the lower leg and calf area.

"Intense Flow" combines the flow of programs A and B. Several chambers always remain inflated at the same time, which prevents reflux. Intense Flow" thus also promotes lactate removal quite intensively. This program is suitable for a fast but intense recovery session.

"Intense x Double x Impulse" combines programs B, C and D, which run sequentially. This gives the user a particularly varied Reboots unit. This can be very pleasant, as the individual leg sections do not experience the pressure as intensively as with one and the same sequence repeated over and over again. An extended regeneration session can be held with this program.

PROGRAM OVERVIEW

H Flow x Double X Intense Flow x Intense



The "Flow x Double x Intense Flow x Intense" program steadily increases the intensity by allowing a lower reflux with each sequence. This means that it is not consistently as intense as program B "Intense", but it is also highly effective and therefore offers beginners an effective alternative. Due to the variety, you can easily hold an entire Reboots unit with program H.

PROGRAM SEQUENCES

#1 RECOVERY IS EVERYTHING

Did you train really hard today? Did you have a big match or race? Then get into your Reboots to boost your recovery!

DURATION: 40 MIN PURPOSE: RECOVERY

10 MIN.	Programm A // FLOW	light pressure
20 MIN.	Programm H // FLOW - DOUBLE - INTENSE FLOW - INTENSE	light - medium pressure
10 MIN.	Programm C // DOUBLE	light pressure

#2 RECOVERY QUICKIE

Do you have a big competition tomorrow? Then quickly get your legs massaged to reach your goals for tomorrow!

DURATION: 20 MIN PURPOSE: QUICK RECOVERY

20 MIN.	Programm B // INTENSE	light - medium pressure
---------	-----------------------	-------------------------

#3 RELEASE YOUR MUSCLES

Treat your legs with some wellness and loosen up your muscles!

DURATION: 40 MIN PURPOSE: MUSCLE RELAXATION

10 MIN.	Programm E // INTENSE FLOW	light pressure
20 MIN.	Programm A // FLOW	light - medium pressure
10 MIN.	Programm F // INTENSE DOUBLE	light pressure

#4 ALL-IN-ONE

Not sure what your legs need the most today? Then we have the right thing for you: A little bit of everything!

DURATION: 40 MIN PURPOSE: REGENERATION & LOOSENING

5 MIN.	Programm E // INTENSE FLOW	light pressure
20 MIN.	Programm G // INTENSE - DOUBLE - IMPULSE	light - medium pressure
10 MIN.	Programm A // FLOW	medium - high pressure
5 MIN.	Programm C // DOUBLE	light pressure

#5 DEEP RELAX

Over and out.

Close your eyes, sit back and relax.

DURATION: 40 MIN PURPOSE: RELAXATION

10 MIN.	Programm A // FLOW	light pressure
20 MIN.	Programm F // INTENSE DOUBLE	light - medium pressure
10 MIN.	Programm B // INTENSE	light pressure

#6 RELAXATION QUICKIE

You need a break from the stressful everyday life? Then just put your feet up for a while and take some time to switch off.

DURATION: 10 MIN PURPOSE: RELAXATION

10 MIN.	Programm H // FLOW - DOUBLE - INTENSE FLOW - INTENSE	medium pressure
----------------	---	-----------------



RECOVERY IS EVERYTHING.